

# The Laminar Wash™ HT2000 System

## Centrifuge-less Sample Preparation for Flow Cytometry and Single Cell Sequencing

The Laminar Wash HT2000 System employs the only suspension-cell sample preparation method that eliminates the centrifuge and the problems it introduces. It is designed to produce the most quantitative and reproducible results for single cell sequencing and flow and mass cytometry.

**Upgrade from the Centrifuge.**

- **Minimize aerosolization** – Reduce the amount of aerosolization compared to a centrifuge and contain infectious samples
- **Increased Cell Retention for low numbers of cells** – For splenocytes and TILs or with rare populations of cells. Reliably high cell retention even with 100's of cells per well
- **Higher viability** – For higher percentages of viable cells after processing compared to centrifugation.
- **Rapid Time to Process** – The system completely washes 96 samples in 4 minutes
- **Better Sequencing Data** – Reduced background with more thorough wash.
- **No Pelleting of Cells** – Reduces doublets, clumping, and clogging.
- **Standardized Results** – Reduces manual pipetting errors and errors associated with multiple personnel changes and locations
- **Higher Stain Index** – For better resolution of populations
- **Cleaner Data** – Improved cell segregation and resolution; Reduces debris and aggregation of cells
- **Small enough to fit inside a biosafety cabinet**


“Laminar Wash technology promotes cell viability and stability and has proven consistency at low cell counts producing reliable and reproducible results.”

*Dr Blanca Ponce-Ngo, Montefiore Medical Center*

Touch-screen interface controls starting volume, flow rate and number of washes.




Laminar Wash Plate with 96 wells



New firmware and PCB allows easy integration with robotic platforms

Optional Buffer Exchanger (BEX) allows user to automatically select up to 8 different solutions from touch-screen interface.


## Laminar Wash Technology Enables Accurate Identification of Tumor Infiltrating Lymphocyte (TIL) Population



(a) Centrifuge wash method – cell loss and mechanical stress through sequential pelleting and resuspension

(b) Laminar Wash method – less tumor debris, but higher retention of TILs and better resolution of populations

**BioPharm** Charles River Laboratories, which uses Laminar Wash HT System by Curiox for TILs



Dimensions	310 mm H x 262 mm W x 302 mm D		
	12.2 in. H x 10.3 in. W x 11.9 in. D		
Voltage Requirements	100 - 240 V	Weight	14 kg
Power Consumption	24V/3.0A		

“With the Curiox Laminar Wash [system], we retain more cells with much less data variation between samples than our centrifuge process.”

*Dr Jorgen Adolffsson, Linkoping University*

Product	Part Number	Description
Laminar Wash™ HT2000 Station 96 Version2	DC-2000-96-01	High throughput 96-well format washing station with GUI
Laminar Wash™ 96-well plate (with a regular lid)	96-DC-CL-05	96-well, coating for flow cytometry assays, sterile
Laminar Wash™ 96 well large volume adaptor	DC-GR-96-05	96 well plate grid to accommodate a larger volume
Buffer Exchanger 5-channel or 10-channel	DC-BX-01-05(10)	Buffer inlet system pedestals to support automation

Go centrifuge-free and accelerate your biology at [curiox.com](https://www.curiox.com)



**CURIOX BIOSYSTEMS INC.**  
400 W Cummings Park  
Suite 4350 Woburn, MA 01801

T 1 781-606-9234  
F 1 650 590 5406  
E [sales@curiox.com](mailto:sales@curiox.com)

Worldwide Contact Info at:  
[www.curiox.com/contact](https://www.curiox.com/contact)  
**Curiox.com**