


# Real-time image-enabled cell sorting

using BD FACSDiscover™ S8 Cell Sorter with BD CellView™ Image Technology

Real-time imaging enhances fluorescence activated cell sorting (FACS) with live visual inspection of target cells and novel gating strategies based on image features. It enables access to detailed information about cells that was previously invisible in traditional flow cytometry experiments and can be used to answer complex biological questions, such as how cells grow, function and interact, or to study exact locations of viruses or proteins within a cell, all at a highly accelerated pace.


## Image Feature Extraction

Image feature extraction with BD CellView™ Image Technology occurs in real-time as raw images are analyzed by firmware components within the cytometer. Image features are numerical measurements of different aspects of the image such as size and shape and are calculated based on the Region of Analysis (ROA). Much like traditional fluorescence, image features can be used to generate dot plots to identify populations of interest which can then be gated and sorted.



Learn how image features are generated using BD CellView™ Image Technology.

Click on the video thumbnail at the right or scan the QR code.



## Image Features

FEATURE AND APPLICATIONS	EXPLANATION	EXAMPLE DIAGRAMS AND HISTOGRAMS
<b>Eccentricity</b> (Scatter and fluorescence) <ul style="list-style-type: none"><li>» Morphology</li><li>» Cluster identification</li><li>» Doublet discrimination</li></ul>	$\text{Eccentricity} = \frac{\text{Long axis}}{\text{Short axis}}$	<ul style="list-style-type: none"><li>■ Low Eccentricity</li><li>■ Int. Eccentricity</li><li>■ High Eccentricity</li></ul>
<b>Radial Moment</b> (Scatter and fluorescence) <ul style="list-style-type: none"><li>» Morphology</li><li>» Cluster identification</li><li>» Doublet discrimination</li></ul>	$\text{Radial Moment} = \text{Median distance of all pixels from the centroid}$	<ul style="list-style-type: none"><li>■ FITC Signal Centrally Concentrated</li><li>■ FITC Signal Radially Spread</li></ul>
<b>Size</b> (Scatter and fluorescence) <ul style="list-style-type: none"><li>» Morphology</li><li>» Label-free sorting</li><li>» Cell cycle</li></ul>	$\text{Size} = \text{Pixel count}$	<ul style="list-style-type: none"><li>■ Small Cells</li><li>■ Large Cells</li></ul>
<b>Total Intensity</b> (Scatter and fluorescence) <ul style="list-style-type: none"><li>» Morphology</li><li>» Imaging FMOs</li></ul>	$\text{Total Intensity} = \text{Additive intensity of all pixels}$	<ul style="list-style-type: none"><li>■ Low Total Intensity</li><li>■ Int. Total Intensity</li><li>■ High Total Intensity</li></ul>

## Image Features (continued)

FEATURE AND APPLICATIONS	EXPLANATION	EXAMPLE DIAGRAMS AND HISTOGRAM
<b>Correlation</b> (Fluorescence: Any two channels) <ul style="list-style-type: none"><li>» Nuclear translocation</li><li>» Co-localization</li><li>» Intracellular trafficking</li></ul>	$\text{Correlation} = \text{Score based on the percentage of pixels from two fluorescence channels that occupy the same space}$	<ul style="list-style-type: none"><li>■ Uncorrelated</li><li>■ Correlated</li></ul>
<b>Delta Center of Mass</b> (Fluorescence: Any two channels) <ul style="list-style-type: none"><li>» Intracellular trafficking</li><li>» Cell-cell interaction</li></ul>	$\text{Delta Center of Mass} = \text{Distance between the centroids of two fluorescence channels}$	<ul style="list-style-type: none"><li>■ Low Delta Center of Mass</li><li>■ Int. Delta Center of Mass</li><li>■ High Delta Center of Mass</li></ul>
<b>Max Intensity</b> (Scatter and fluorescence) <ul style="list-style-type: none"><li>» Punctate fluorescence</li><li>» Phagocytosis assay</li><li>» Cell cycle analysis</li></ul>	$\text{Max Intensity} = \text{Brightest pixel}$	<ul style="list-style-type: none"><li>■ Dim Fluorescent Cells</li><li>■ Int. Fluorescent Cells</li><li>■ Bright Fluorescent Cells</li></ul>
<b>Diffusivity</b> (Scatter and fluorescence) <ul style="list-style-type: none"><li>» Punctate fluorescence</li><li>» Cell morphology</li><li>» Phagocytosis</li></ul>	$\text{Diffusivity} = \frac{\text{Fluorescence-A}}{\text{Max Intensity}}$	<ul style="list-style-type: none"><li>■ Low Diffusivity</li><li>■ Int. Diffusivity</li><li>■ High Diffusivity</li></ul>

BD flow cytometers are Class 1 Laser Products.  
For Research Use Only. Not for use in diagnostic or therapeutic procedures.  
Distributed by BD Life Sciences, 155 North McCarthy Blvd, Milpitas, CA 95035, U.S. | bdbiosciences.com  
BD, the BD Logo, BD CellView and BD FACSDiscover are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2024 BD. All rights reserved. BD-136898 (v1.0) 0924

