



Objective.

Protocol for cell lysis, designed for whole proteome characterization by LC-MS/MS. The protocol uses nonionic detergent (NP-40) and is not suitable if your goal is to specifically analyze membrane proteins.

Please follow this protocol and provide the Proteomics Resource Center with 10-50ug of quantified protein.

The protocol is provided by Maria Passarelli.

Lysis buffer	Concentration	
Tris-HCl	20mM	
KCl	100mM	
NP-40	0.5%	Do <u>not</u> substitute with SDS
cOmplete Protease Inhibitor (EDTA free)	1 pellet per 10mL	

Other equipment

Centrifuge

Sonicator

BCA assay

Equipment for running SDS-PAGE



Procedure – Cell lysis

Examples, Tricks & Comments

-
- 1 Prepare lysis buffer as described above

 - 2 Resuspend cell pellet in appropriate volume of lysis buffer (400 uL minimum for sonication) in 1.5 mL Eppendorf tube

 - 3 Vortex for 5 seconds

 - 4 Sonicate samples in ice (amplitude: 40%) at 5 second pulses with 30 second breaks, 4 pulses total

 - 5 Centrifuge samples at maximum speed for 10 minutes at 4C.

 - 6 Transfer supernatant to new tubes.

Procedure – Quality Control

-
- 1 Withdraw an aliquot from all samples and measure concentration using BCA or comparable assay.

 - 2 Withdraw 5 ug from each sample and run SDS-PAGE. Stain with Colloidal blue. If all lanes look similar, your sample is ready to submit. *If lanes do not look similar, it may be caused by contamination of a single high-abundant protein. Repeat lysis.*
-