

Guidelines for the Hands-on Training on BD AccuriC6

1. **Ice bucket** with samples covered with aluminum foil
Note: All the samples (unless otherwise specified by the protocol) should remain at +4°C protected from light
2. Plenty of the **unstained cells**. Cell concentration to be $1-5 \times 10^6/\text{ml}$. 1 ml volume should be enough. Useful Tips to Improve Sample Quality could be found at <http://www.rockefeller.edu/fcrc/tips/sampleprep>
Note: These unstained cells should fully match the sample origin and preparation to be used in the future data acquisition on Accuri C6.
3. **Cell samples stained with different concentrations of antibody** in order to find the proper titer for each antibody (if applicable):
 - a. For the antibody previously not used in the lab – bring "wide" titration.
For example: 3x; x; 1/3x; 1/9x; 1/27x; 1/81x (where "x" is the concentration suggested by the vendor)
 - b. For the antibody previously used in the lab – bring "narrow" titrations.
For example: 3y; y; 1/3y; 1/9y (where "y" is the concentration suggested by the Labmates)
 - c. For Live Dead Fixable Dyes - bring "narrow" titration.
For example: 2z; z; 1/2z (where "z" is the concentration suggested by the vendor). Make sure to follow the vendor's staining protocol (time, temperature, special buffer)
Note: Cell concentrations to be $1-5 \times 10^6/\text{ml}$. 500 μl volume should be enough.
4. 15 ml tube with the **"FACS Buffer"** (<http://www.rockefeller.edu/fcrc/tips/cellsorting>) to use in case we need to dilute the samples or titer non-fixable Dead Cell Exclusion Dye (DCE).
5. 1.5 ml of the 1 $\mu\text{g}/\text{ml}$ **stock of the non-fixable DCE** (for example 7AAD, PI, ToPro3, etc.) to prepare titration and stain samples.
Note: These dyes remain in the "FACS Buffer" with single cell suspension (no wash step).
6. **Automatic pipets** (20, 200 and 1000 μl).
Note: FCRC doesn't lend pipets in the Analysis room, but offers tips.
7. **Sharpie marker** (to label tubes).