



















Aurora R0018





Main Fluorochrome Signatures
@ Cytek Assay Settings

Violet Dyes





Red Dyes



Peak Channel for Tested Fluorochromes

Fluorochrome	Primary Channel	Fluorochrome	Primary Channel	Fluorochrome	Primary Channel
BV421	V1	BB515	B1	APC	R1
Super Bright 436	V2	FITC	B2	APC-Cy5.5	R3
eFluor450	V3	Alexa532	B3	Alexa Fluor700	R4
BV510	V6	PE	B4	APC-H7	R8
BV570	V8	PE-CF594	B6		
BV605	V10	PE Cy5	B8		
BV650	V11	PE-Cy5.5	B9		
BV711	V13	PerCP-eFluor710	B10		
BV785	V15	PE Cy7	B13		

