

# Fluorochrome Peak Channels

UV Excited Fluors	Peak Channel
BUV395	UV2
BUV496	UV7
BUV563	UV9
BUV661	UV11
BUV737	UV14
BUV805	UV16

Violet Excited Fluors	Peak Channel
BV421	V1
Alexa 405, SuperBright 436	V2
eFluor 450, VioBlue, Pacific Blue	V3
BV480	V4
eFluor 506	V5
BV510, VioGreen	V7
BV570, Pacific Orange	V8
BV605, SuperBright 600, Qdot 605	V10
BV650, SuperBright 645, Qdot 655	V11
BV711, SuperBright 702, Qdot 705	V13
BV750	V14
BV785, BV786, Qdot 800	V15

Blue Excited Fluors	Peak Channel
Vio 515, sVio 515, BB515	B1
Alexa Fluor 488, FITC, VioBright FITC	B2
Alexa Fluor 532	B3
PerCP	B8
PerCP-Cy5.5, BB700	B9
PerCP-eFluor 710, PerCP-Vio 700	B10

Yellow Green Excited Fluors	Peak Channel
PE	YG1
PE-Dazzle 594, PE-CF 594, PE-TexasRed, PE-eFluor 610, PE-Alexa Fluor 610	YG3
PE-Cy5	YG5
PE-Cy5.5, PE-Alexa Fluor 700	YG7
PE-Cy7, PE-Vio 770	YG9

Red Excited Fluors	Peak Channel
APC	R1
Alexa 647, Vio 667, sVio 667, eFluor 660	R2
Alexa 700, APC-R700	R4
APC-Alexa 750, APC/Fire 750, APC-Cy7, APC-Vio 770, APC-eFluor 780, APC-H7	R7

